

Hemostasis test results, particularly global assays such as prothrombin time (PT) or activated partial thromboplastin time (aPTT), are known to vary with the reagents and analyzers used.⁷ The Perinatal and Pediatric Hemostasis Subcommittee of the SSC of the ISTH recommends that each laboratory define age-dependent reference ranges using their own methodologies.⁸

To assist labs in this endeavor, IL collaborated with six hospital laboratories to define reference ranges for the use of HemosIL assays in children.³ The results from this collaboration³ confirmed that most Hemostasis test results are highly dependent on age, and that age-specific reference ranges must be used to ensure proper evaluation.

References

1. Appel IM, *et al.* Age dependency of coagulation parameters during childhood and puberty. *J Thromb Haemost.* 2012;10:2254–63.
2. Andrew M, *et al.* Developmental Hemostasis: Relevance to Hemostatic Problems During Childhood. *Semin Thromb Hemost.* 1995;21:341–56.
3. Toulon P, *et al.* Age dependency for coagulation parameters in pediatric populations: results of a multicenter study aimed at defining the age-specific reference ranges. *Thromb Haemost.* 2016 Mar 17;116(1) [Epub ahead of print].
4. Toulon P, *et al.* Age dependency of coagulation parameters during childhood. Poster presented at 57th ASH Annual Meeting, Orlando, FL, December 2015.
5. Trimm A, *et al.* Comparison of automated von Willebrand factor activity assays. *Thrombosis Research.* 2015;135:684–91.
6. Costa-Pinto J, *et al.* Diagnosis of inherited von Willebrand disease: comparison of two methodologies and analysis of the discrepancies. *Haemophilia.* 2014;20(4):559–67.
7. Lippi G, *et al.* Coagulation testing in pediatric patients: the young are not just miniature adults. *Semin Thromb Haemost.* 2007;33:816–20.
8. Ignjatovic V, *et al.* On behalf of the Perinatal and Paediatric Haemostasis Subcommittee of the SSC of the ISTH. *J Thromb Haemost.* 2012;10:298–300.

► HEMOSTASIS INNOVATION IS HERE.

The Instrumentation Laboratory logo, ACL, ACL AcuStar, ACL Elite, ACL TOP, HemoCell, HemoHub, HemosIL, ReadiPlasTin, RecombiPlasTin, SynthAFax and SynthASil are trademarks of Instrumentation Laboratory Company and/or one of its subsidiaries or parent companies and may be registered in the United States Patent and Trademark Office and in other jurisdictions. All other product names, company names, marks, logos and symbols are trademarks of their respective owners.

©2017 Instrumentation Laboratory. All rights reserved.

For more information, contact your local Werfen sales representative or distributor.

CODE

Optimizing Pediatric Hemostasis Testing.

About this guide

Diagnostic testing is critical to ensuring optimal patient management of hemorrhagic and thrombotic diseases in children.^{1,2} This guide presents the reference ranges of diagnostic assays for pediatric patients (ages 15 days to 17 years), as determined by a multi-center study using ACL TOP® Hemostasis Testing Systems.³

Age dependency of Hemostasis parameters during childhood³⁻⁶

Pediatric reference ranges for HemosIL Assays on ACL TOP Hemostasis Testing Systems

Test	HemosIL Assay	15 Days-4 Weeks	1-5 Months	6-11 Months	1-5 Years	6-10 Years	11-17 Years
PT (sec)	RecombiPlasTin® 2G	11.2 [9.5 - 12.6]	11.0 [9.7 - 12.8]	11.0 [9.8 - 13.0]	11.3 [9.9 - 13.4]	11.7 [10.0 - 14.6]	11.8 [10.0 - 14.1]
aPTT (sec)	SynthASil®	35.4 [27.6 - 45.6]	33.5 [24.8 - 40.7]	32.4 [25.1 - 40.7]	31.6 [24.0 - 39.2]	31.6 [26.9 - 38.7]	31.0 [24.6 - 38.4]
aPTT (sec)	aPTT-SP	39.0 [33.2 - 45.6]	33.3 [25.0 - 43.3]	33.3 [25.0 - 43.3]	32.4 [25.7 - 38.4]	32.8 [25.5 - 42.4]	32.6 [26.1 - 47.4]
Fibrinogen (mg/dL)	Fibrinogen-C	254 [143 - 402]	226 [150 - 376]	233 [157 - 360]	273 [188 - 413]	278 [189 - 475]	266 [177 - 420]
Fibrinogen (mg/dL)	Q.F.A. Thrombin	240 [136 - 300]	210 [141 - 437]	230 [148 - 367]	260 [164 - 497]	276 [171 - 537]	248 [168 - 529]
FII (%)**	Factor II Deficient Plasma	56.3 [44.8 - 74.3]	75.0 [46.7 - 110.6]	91.5 [73.9 - 117.2]	99.0 [49.4 - 130.0]	90.0 [68.4 - 132.0]	93.5 [47.6 - 119.2]
FV (%)**	Factor V Deficient Plasma	100.0 [69.0 - 123.7]	99.5 [59.5 - 147.0]	102.0 [59.0 - 159.8]	110.5 [73.2 - 188.1]	101.0 [82.0 - 140.6]	97.0 [61.7 - 124.8]
FVII (%)**	Factor VII Deficient Plasma	75.6 [55.0 - 108.0]	88.0 [43.0 - 141.1]	88.0 [55.2 - 128.0]	82.0 [47.8 - 124.2]	77.0 [55.0 - 135.4]	81.5 [55.4 - 133.1]
FVIII (%)*	Factor VIII Deficient Plasma	95.5 [65.2 - 153.4]	84.5 [50.3 - 187.3]	75.0 [53.4 - 132.2]	95.0 [59.0 - 167.0]	86.5 [60.6 - 154.4]	93.0 [42.8 - 154.6]
FIX (%)*	Factor IX Deficient Plasma	43.5 [30.0 - 77.0]	53.0 [29.0 - 105.1]	76.5 [50.5 - 106.8]	84.0 [52.6 - 128.9]	80.0 [55.3 - 156.0]	96.5 [60.2 - 138.4]
FX (%)**	Factor X Deficient Plasma	85.0 [66.0 - 92.0]	89.0 [67.5 - 122.2]	100.0 [75.8 - 134.4]	99.0 [59.7 - 152.8]	99.0 [71.3 - 161.5]	93.0 [64.0 - 130.5]
FXI (%)*	Factor XI Deficient Plasma	56.0 [32.9 - 75.0]	64.0 [27.6 - 126.4]	86.0 [60.9 - 125.6]	92.0 [58.0 - 154.0]	83.0 [31.8 - 154.0]	84.0 [55.4 - 139.4]
FXII (%)*	Factor XII Deficient Plasma	69.2 [25.0 - 81.0]	76.0 [38.0 - 136.6]	109.2 [48.0 - 156.1]	107.0 [50.0 - 174.7]	83.7 [49.4 - 153.5]	91.7 [49.4 - 153.5]
FXIII (%)	Factor XIII Antigen	86.0 [78.4 - 100.0]	82.9 [55.3 - 133.2]	92.0 [51.1 - 136.8]	97.4 [49.4 - 137.2]	96.5 [53.5 - 142.4]	106.0 [64.4 - 133.1]
vWF:GP1bR (%)	von Willebrand Factor Ristocetin Cofactor Activity	99.6 [87.8 - 121.5]	89.0 [33.2 - 154.1]	67.1 [37.1 - 118.6]	83.3 [40.8 - 131.8]	89.1 [42.1 - 162.6]	92.8 [45.0 - 139.1]
vWF:Ab (%)	von Willebrand Factor Activity	121.5 [73.7 - 188.9]	104.0 [40.9 - 191.0]	86.0 [42.7 - 176.0]	82.4 [43.6 - 155.8]	83.0 [41.2 - 128.9]	83.5 [54.0 - 136.9]
vWF:Ag (%)	von Willebrand Factor Antigen	163.3 [46.0 - 219.5]	101.5 [35.5 - 217.0]	78.6 [48.4 - 199.4]	89.1 [41.0 - 185.7]	80.0 [44.8 - 141.1]	92.0 [55.6 - 123.4]
Antithrombin (%)	Liquid Antithrombin	41.0 [32.8 - 62.8]	80.1 [29.0 - 120.0]	96.0 [63.0 - 121.8]	96.5 [60.5 - 128.3]	97.0 [64.2 - 136.4]	97.0 [69.1 - 135.9]
PC chromogenic (%)	Protein C	39.1 [27.2 - 48.0]	51.2 [22.8 - 95.0]	79.9 [46.6 - 150.9]	92.6 [59.1 - 147.5]	100.5 [45.9 - 153.5]	99.0 [72.3 - 155.1]
PC clotting (%)	ProClot	37.5 [29.7 - 114.6]	82.0 [28.1 - 127.8]	85.0 [43.7 - 151.3]	86.3 [61.0 - 143.5]	91.0 [39.3 - 170.3]	95.1 [65.8 - 126.6]
PS Free Ag (%)	Free Protein S (antigenic immunoassay)	83.8 [61.0 - 108.0]	94.9 [48.0 - 126.5]	86.0 [63.0 - 138.9]	86.4 [53.0 - 134.9]	95.1 [61.5 - 141.7]	93.5 [61.4 - 130.7]
PS clotting (%)	Protein S Activity	90.1 [29.0 - 115.2]	81.6 [33.3 - 153.9]	88.3 [51.8 - 138.4]	97.6 [60.2 - 148.8]	104.8 [66.5 - 161.5]	99.3 [52.5 - 147.1]
Plasminogen (%)	Plasminogen	52.6 [41.0 - 82.7]	69.2 [37.6 - 109.6]	80.7 [49.3 - 126.4]	91.8 [59.6 - 178.0]	92.0 [52.4 - 158.1]	91.8 [58.1 - 130.6]
D-dimer (ng/mL FEU)	D-Dimer HS 500	530 [445 - 1200]	515 [90 - 878]	270 [133 - 844]	280 [88 - 780]	275 [60 - 567]	245 [69 - 580]

Pediatric claims have not been reviewed by regulatory agencies. Listed products should be used per labeled claims. Units of measure from publication converted to those listed on Product Insert Sheets.

*Used with SynthASil.

**Used with RecombiPlasTin 2G.